

Characteristics of the New Gamma Generator:

Comfort:

1. Small and handy
2. Intuitive handling and quick reference for quicker understanding
3. Compatible with other electrodes available from us
4. LCD display with icons
5. Complete accessories

Technology:

6. Generates sine or square wave
7. Can be used as a Zapper
8. Generates a frequency between 1Hz and 1MHz, which covers all Clark and Rife frequencies
9. Exact signal
10. Programs with up to 16 program steps
11. Programs can be stored on external storage (Driver) (see table on reverse side of this info sheet)
12. Offset adjustable (sine wave generator)
13. Contact check can be deactivated for use without load (e.g. with Zincrometer, Zappicator)
14. TENS Pain Treatment with bioTENSion unit (US: pending FDA approval)

Design:

15. Attractive design
16. Practical carrying case

Compatibility:

15. Attach to Zincrometer
16. Attach to plate zapping plates
17. Connect to Zappicator
18. Produce homeographic copies
19. Run Zapper Drivers
20. Run frequency programs from BioStream
21. Store frequency programs for Gamma or BioStream

Accessories:

- 1 empty sine and 1 empty square Driver
- 1 pair of gel electrodes
- 1 pair of wrist bands
- 1 electrode cable
- 1 cable with alligator clip for connection to Zincrometer, plate zapping plates, plates for homeographic copies
- 2 AA batteries

